




technoclone



FACTOR DEFICIENT PLASMA


Factor II, V, VII, VIII, IX, X, XI, XII,
Fitzgerald Trait and Fletcher Trait deficient plasma

Factor VIII, IX, X, XI, XII deficient plasma is available native and immunodepleted!

A clear view into the future!

FACTOR ASSAYS

Factors of the intrinsic Pathway are tested in one-stage clotting assays with activated Partial Thromboplastin Time (aPTT) reagents using factor VIII, IX, XI and XII deficient plasma as well as Fitzgerald Trait and Fletcher Trait plasma.



Factors of the extrinsic and Common Pathway are tested in one-stage clotting assays with Prothrombin Time (PT) reagents using plasma deficient of Factor II, V, VII and X.

Product	Description	REF	Package
Factor II deficient plasma	Lyophilized Factor II deficient plasma, immunodepleted	5114008	5 x 1 ml
Factor V deficient plasma	Lyophilized Factor V deficient plasma, immunodepleted	5134004	5 x 1 ml
Factor VII deficient plasma	Lyophilized Factor VII deficient plasma, immunodepleted	5144015	5 x 1 ml
Factor VIII deficient plasma	Lyophilized Factor VIII deficient plasma, immunodepleted	5154002	5 x 1 ml
		5154004	50 x 1 ml
	Lyophilized Factor VIII deficient plasma, native	5154007	5 x 1 ml
		5154016	50 x 1 ml
Factor IX deficient plasma	Lyophilized Factor IX deficient plasma, immunodepleted	5164003	5 x 1 ml
		5164004	50 x 1 ml
	Lyophilized Factor IX deficient plasma, native	5164008	5 x 1 ml
		5164016	50 x 1 ml
Factor X deficient plasma	Lyophilized Factor X deficient plasma, immunodepleted	5174006	5 x 1 ml
	Lyophilized Factor X deficient plasma, native	5174004	5 x 1 ml
Factor XI deficient plasma	Lyophilized Factor XI deficient plasma, native	5184004	5 x 1 ml
	Lyophilized Factor XI deficient plasma, immunodepleted	5184006	5 x 1 ml
Factor XII deficient plasma	Lyophilized Factor XII deficient plasma, native	5194008	5 x 1 ml
	Lyophilized Factor XII deficient plasma, immunodepleted	5194007	5 x 1 ml
Fitzgerald Trait plasma	Lyophilized high molecular weight kininogen deficient plasma, immunodepleted	5204006	2 x 1 ml
Fletcher Trait plasma	Lyophilized prekallikrein deficient plasma, native	5205006	2 x 1 ml

