

technoclone

CEVERON[®] MFU 500

Microparticle filtration unit



A clear view into the future!

Ceveron[®] MFU 500


Ceveron[®] MFU 500 enables easy, quick and standardised removal of *MICROPARTICLES* from platelet poor plasma.




EASY & QUICK



Pipette sample into the filter plate



Push the "START" button and filtration starts automatically



After filtration is terminated push the "VENTILATION" button



Filtered sample can be collected from the receiver plate

SMALL & COMPACT




Filtration is as effective as ultra-centrifugation to remove thrombogenic MP's from PPP.

Dimensions MFU 500:

Length: 26.00 cm
Width: 20.00 cm
Height: 16.05 cm
Weight: 3.58 kg

STANDARDIZED



Product	Description	REF	Package
Ceveron [®] MFU 500	Microparticle filtration unit	9821010	1 instrument

